

Low Intrinsic Clearance determination using Primary Human Hepatocytes in monoculture compared to co-culture with non-parenchymal stromal cells.

Petter Svanberg¹, Britta Bonn¹, Malin Darnell¹, Ia Hultman¹ and Ken Grime¹

¹AstraZeneca R&D Mölndal, Sweden



petter.svanberg@astrazeneca.com

Abstract

- Definition of intrinsic clearance (CL_{int}) values of less than $1\mu\text{L}/\text{min}/\text{million}$ human hepatocytes is challenging since this implies a drug half-life of 12 hours or less in an hepatocyte CL_{int} assay
- Drug metabolizing enzyme (DME) activities decline to low levels after 4 hours in primary human hepatocyte suspension incubations¹ but pooled multi-donor platable hepatocytes in monoculture on collagen-coated plates have shown slower decline of activity for important DMEs and stability is maintained for 8-12 hours of incubation²
- H μ REL Corporations novel co-culture of hepatocytes with stromal cells has shown sustained activity of DMEs up to 2-3 weeks in culture³
- We compared the two formats of cultured primary human hepatocytes (5-donor pool, XenoTech) & the potential to provide improved clearance prediction for drugs with high metabolic stability
- Enzyme activity changes (EAC) during the incubation were monitored by determining rate of metabolite formation for DME selective probe substrates at various times through the incubations
- Results demonstrate that hepatocytes plated in both monoculture (10 hour incubation) and in co-culture (72 hour incubation) can robustly produce CL_{int} values less than $1\mu\text{L}/\text{min}/\text{million}$ cells
- Both systems predict human clearance well but the combination of high hepatocyte density & DME activity stability up to 72 hours, as well as assay-ready plate delivery, makes the H μ REL co-culture system a preferable method for determination of low CL_{int}

Objectives

- To compare low CL_{int} determination & prediction of human *in vivo* clearance using pooled platable hepatocytes in monoculture & co-culture of hepatocyte/non-parenchymal stromal cells.

Methods

- CL_{int} was determined for a set of substrates (Table 1) with primary human hepatocytes in monoculture for 10 hours as well as in co-culture for 72 hours. EAC was monitored with a CYP-selective substrate cocktail (Midazolam, Diclofenac, Bufuralol, S-Mephenytoin, Phenacetin) assessing formation of DME specific metabolites for 1 h at different time points after plating.

Monoculture

Pooled human hepatocytes (lot 1310168) were purchased from XenoTech, thawed and seeded into collagen-I-coated 24 well plates for 4 hours, according to XenoTechs plating protocol. Incubations were then started.

Co-culture

Pooled human hepatocytes (lot 1310168) were co-cultured with non-parenchymal stromal cells for six days at H μ REL, US before subsequent shipping of assay ready 24 well plates to AstraZeneca labs in Sweden. Plates were acclimatized for 24 hours before start of incubation at day 7 after seeding.

- Fresh culture media containing typically $1\mu\text{M}$ substrate was added to cells & aliquots of media were withdrawn from each well during incubations in order to determine CL_{int} .
- Substrates & metabolites were quantified using LCMSMS. CL_{int} was calculated from 'parent drug-loss' plots & scaled into *in vivo* CL_{int} .
- Mathematical scaling was performed using the following parameters & method: 120 million hepatocytes per g liver, 24 g liver per kg bodyweight. Scaled *in vivo* CL_{int} values were corrected with a regression line off-set approach to predict *in vivo* hepatic CL_{int} ⁴. The scaled *in vivo* CL_{int} values were compared to back-scaled (well-stirred model) *in vivo* CL_{int} from human plasma clearance.

Results and Discussion

- EAC (decline) for major CYPs gives a window of approximately 10 hours for log-linear compound disappearance (metabolism) from monocultures. In contrast, EAC (decline) in the co-culture system is considerably more stable over 72 hours (data not shown).
- CL_{int} values in the range of 0.2 to $34\mu\text{L}/\text{min}/\text{million}$ hepatocytes were determined. A suspension assay format gave calculable CL_{int} values for only 11 of the 21 compounds whereas 17/21 values were determinable from the monoculture data and 19/21 from the co-culture assay (Figure 1, Table 1). Typical compound depletion profile used for CL_{int} determinations are shown in Figure 2.
- Monoculture and co-culture predicted *in vivo* CL_{int} to within 2 to 3-fold for the majority of tested compounds: Monoculture average fold difference (AFD) of 1.8; Co-culture AFD of 2.6 (Table 1).

Figure 1.

Success rate in determination of CL_{int} in different hepatocytes CL_{int} assays.





Figure 2.

Disopyramide disappearance from media in H μ RELs co-culture.



Conclusions

- Plated human hepatocytes in 10 hour monoculture incubations & in stromal cell co-culture for 72 hours predict *in vivo* clearance to within 2 to 3-fold for compounds metabolised by key DMEs
- Both systems show high success rate for determination of CL_{int} in the 0.2 to $34\mu\text{L}/\text{min}/\text{million}$ hepatocyte range (CL_{int} below 1 is readily measurable)
- H μ RELs assay-ready plates improve success rate for CL_{int} determination

Table 1. Compound information, observed CL_{int} -values and *in vivo*-prediction data. (Between 2 and 6 replicate incubations with *in vivo* predicted compounds).

Test compound properties			Observed <i>in vitro</i> CL_{int}			Observed <i>in vivo</i> CL_{int}	Predicted <i>in vivo</i> CL_{int}			
Compound	Ion class	DME	Co-culture	Monoculture $\mu\text{L}/\text{min}/\text{million}$ cells	Suspension	$\text{mL}/\text{min}/\text{kg}$	Co-culture		Monoculture	
							$\text{mL}/\text{min}/\text{kg}$	Fold ^a	$\text{mL}/\text{min}/\text{kg}$	Fold ^a
AZ001	Neutral		n.v.	0.4	n.v.	2.1 ^b	n.v.	n.v.	1.9	-1.1
Theophylline	Neutral	1A2	n.v.	n.v.	n.v.	1.4 ⁵				
Disopyramide	Base	3A4	0.4	0.2	1.5	0.9 ⁶	1.5	1.6	0.9	-1.1
S-warfarin	Neutral	2C9>3A4	0.7	n.v.	n.v.	0.1 ⁵	0.4	6.5	n.v.	n.v.
AZ002	Acid		0.8	n.v.	n.v.					
AZ003	Base		0.8	n.v.	n.v.					
AZ004	Neutral		0.8	0.4	n.v.					
Metoprolol	Base	2D6	0.9	1.4	n.v.	37.1 ⁶	8.8	-4.2	12.3	-3.0
Diazepam	Neutral	2C19>3A4	1.3	0.8	n.v.	0.5 ⁵	1.2	2.3	0.8	1.5
AZ005	Base		1.4	1.6	n.v.					
Imipramine	Base	2C9>2D6>3A4>1A2	1.7	8.6	4.5	41.5 ⁵	15.3	-2.7	50.9	1.2
AZ006	Neutral		1.8	1.0	n.v.					
AZ007	Neutral		1.9	n.v.	n.v.					
AZ008	Neutral		2.0	1.3	4.2					
AZ009	Acid	3A4	3.3	1.2	3.6	4.7 ^b	2.3	-2.1	1.1	4.4
Ketoprofen	Acid	UGT2B2	4.3	3.9	6.6	2.4 ⁵	8.1	3.3	7.5	3.1
Verapamil	Base	3A4	4.9	11.1	18.5	39.7 ⁵	17.2	-2.3	31.9	-1.2
Sildenafil	Base	3A4 > 2C9/2C19	6.2	7.0	9.9	8.6 ⁵	7.2	-1.2	7.9	1.1
Ketanserin	Base	non-CYP AKR/CR	14.1	8.3	33.1	10.1 ⁶	15.3	1.5	14.7	1.5
AZ010	Base		27.3	28.4	5.9					
Carvedilol	Base	2D6>2C9	34.2	26.3	37.2	15.4 ⁵	19	1.2	15.6	1.0

^a Fold difference predicted versus observed *in vivo* CL_{int}

^b In-house clinical data

References

1. Smith M. C et al. Journal of Pharmaceutical Sciences, 101 (10), pp 3989-4002, 2012.
2. Svanberg P et al. Poster at MDO symposium, Stuttgart, 2014.
3. <http://hurelcorp.com/products-services/>
4. Sohlenius-Sternbeck, A-K.; et al. Xenobiotica 42(9):841-85, 2012.
5. Halifax D1, Foster JA, Houston JB. Pharm Res. Oct;27(10):2150-61, 2010
6. Obach RS1, Lombardo F, Waters NJ. Drug Metab Dispos. Jul;36(7):1385-405, 2008.