



HμREL  
CORPORATION

# HμREL*flux*<sup>TM</sup>

## Biliary Efflux Assay Kit

---

### Frequently Asked Questions (FAQ's)

---

1) **What is HμREL*flux*<sup>TM</sup>?**

HμREL*flux*<sup>TM</sup> is a novel, patent-pending method for assessing a compound's biliary disposition.

2) **Why is biliary disposition important?**

Transporters can play an important role in the clearance and excretion of the parent drug and/or its metabolite(s) into the bile canaliculus. In-vitro methods for accurately predicting such dispositions are of great need.

3) **What data does the HμREL*flux*<sup>TM</sup> assay return to the user?**

With this assay one can measure:

- Compound clearance from the dosing solution
- Biliary accumulation
- Intracellular accumulation

With these data one can predict a compound's:

- Intrinsic clearance
- Uptake rate
- Biliary clearance
- Biliary excretion index

4) **Does HμREL*flux*<sup>TM</sup> use gel overlays?**

No. Gels can be cumbersome to apply and may cause non-specific binding of your compounds. HμREL<sup>®</sup> primary hepatic co-cultures readily form bile canaliculi, there is no need to use special "Certified" cells

5) **How does the HμREL<sup>®</sup> model achieve the levels of hepatic function and morphology necessary to form bile canaliculi?**

HμREL<sup>®</sup> uses its own patented primary hepatocyte co-cultures which allow the hepatocytes to polarize and form stable and enduring biliary canaliculi networks and their associated transporters for over 30 days.

**6) How does HμREL® measure the relatively small quantities of bile in the cultures?**

In order to accurately quantify the biliary clearance of a compound we:

- a) Take all of our measurements in one well and assess the biliary portion and cellular portion from the same single well.
- b) Take accurate measurement of small quantities by LC/MS.

**7) Do the compounds being assessed need to be radio labeled?**

No, radio or fluorescent labeling is not required. LC/MS can provide accurate measurement without need to label compounds.

**8) Is it possible to generate a negative biliary excretion index (“BEI”)?**

No. Because HμRELflux™ takes direct measurement of both the biliary and, separately, the intracellular portions—all from within a single well of hepatocytes, with no arithmetic comparison of results obtained from multiple wells—it is impossible to generate a negative BEI number.

**9) How does HμREL® account for basolateral transport during the canalicular disruption phase of the experiment?**

While it is certainly possible for a compound to be effluxed via basolateral transporters in vitro; we believe that this parameter is negligible for most compounds and have engineered our assay to minimize the effect that this phenomena may have on your biliary excretion readings.

**10) Is HμRELflux™ available as a kit?**

Yes, The HμRELflux™ biliary efflux assay kit is air-shipped by overnight courier to your lab, where it arrives ready for immediate, “plug-and-play” use after a brief period of a few hours’ acclimation. No plating and only minimal tissue culture is required of the customer.

**11) What does the HμRELflux™ biliary efflux assay kit contain?**

Media (both maintenance and dosing media), HμREL® primary hepatocyte co-cultures, all reagents necessary to perform the HμRELflux™ assay, detailed protocol, unpacking protocol and certificate of analysis will be included in the HμRELflux™ kit.

**12) Is HμRELflux™ also available on a contract research services basis?**

Yes. HμREL® expert technical staff will perform the HμRELflux™ assay on your compound(s) if you prefer to outsource this workflow. Please contact HμREL® for more information on HμRELflux™ contract research services.

**13) Liver models of which species have been characterized for HμRELflux™?**

HμREL® primary hepatic co-cultures are available in a variety of species. Currently, the HμRELhuman™, HμRELprimate™, HμRELDog™, and HμRELrat™ patented models have been qualified using the HμRELflux™ method.

**14) Which uptake transporters have been assessed using the HμRELflux™ method?**

Probe substrates	Uptake transporters
Taurocholic acid (2μM)	NTCP, OATPs
Rosuvastatin (2μM)	OATPs
Digoxin (2μM)	OATP-8
Estradiol-17b-D-glucuronide (2μM)	OATPs
Pravastatin (5 μM)	OATPs

**15) Which efflux transporters have been assessed using the HμRELflux™ method?**

Probe substrates	Uptake transporters
Taurocholic acid (2μM)	BSEP
Rosuvastatin (2μM)	BCRP/MRP2
Digoxin (2μM)	P-gp, mdrp-3
Estradiol-17b-D-glucuronide (2μM)	MRP2
Pravastatin	BCRP/MRP2

**16) How long should you wait after receiving the cells before beginning your experiment?**

You can start your experiments using HμREL®'s primary hepatic co-cultures after a minimum four-hour re-acclimation period in an incubator, after replacement of the shipping media with the supplied maintenance media that is provided as a component of the kit. Alternatively, you can wait as long as three days after receiving your HμREL® co-cultures before beginning your experiments, with no diminution of hepatic function.

**17) How long will an experiment take using the HμRELflux™ method?**

Our customers have conducted their experiments over various time courses. Given the robustness of the HμREL® system, one has the ability to evaluate a variety of different rates, from a minimum of one minute to a maximum of four days. The HμRELflux™ method allows for easy biliary uptake/efflux evaluation, even of compounds with properties that make evaluation difficult when using traditional transporter evaluation methods.

**18) Have in vitro-to-in vivo (“IVIVC”) correlations been made for results obtained using HμRELflux™?**

Yes. Using the HμRELrat™ hepatic co-culture model, in vitro-to-in vivo correlation was within 3-fold for all instances.

**For more information:**

Email the company at [service@hurelcorp.com](mailto:service@hurelcorp.com)

Or call: (732) 253-0295